EBA 20

Ⓧ DE	Bedienungsanleitung.....	4
Ⓧ EN	Operating Instructions.....	14
Ⓧ FR	Mode d'emploi	24
Ⓧ IT	Istruzioni per l'uso	34

**EG-Konformitätserklärung
EC Conformity Declaration
Déclaration de conformité CE
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

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Contents

1	Use according to specification	15
2	Residual risks	15
3	Notes on safety.....	15
4	Symbol meanings	16
5	Delivery checklist.....	16
6	Space requirement	17
7	Connection to the mains	17
8	Commissioning.....	17
9	Opening and closing the lid	17
9.1	Opening the lid	17
9.2	Closing the lid.....	17
10	Arming the rotor.....	17
11	Control panel	18
11.1	Speed area	18
11.2	Rotation area	18
11.3	Time area.....	18
11.4	Key area	18
12	Adjustable parameters.....	19
13	Preselecting centrifuging parameters or changing them during operation	19
13.1	Speed	19
13.2	- of denser substances.....	19
13.3	Run time	19
13.3.1	Continuous operation.....	19
13.3.2	Pulsed operation.....	19
13.4	Brake adjustment.....	20
13.5	Start centrifuging run.....	20
13.6	End centrifuging run.....	20
14	Calculating rotational speed RPM and relative centrifugal force RCF	20
15	Emergency release.....	21
16	Maintenance and servicing.....	21
16.1	Centrifuge	21
16.2	Rotors and Attachments	21
16.3	Autoclaving	22
16.4	Centrifuge containers.....	22
17	Faults.....	22
18	Change mains input fuse.....	23
19	Acceptance of the centrifuges for repair	23
20	Anhang / Appendix	44
20.1	Technische Daten / Technical specification	44
20.2	Rotoren und Zubehör / Rotors and accessories.....	45

1 Use according to specification

The machine presented here is a medical product (laboratory centrifuge) according to the IVD guideline 98/79/EG. The centrifuge is used to separate substances or substance mixtures with a density of max. 1.2 kg/dm³. This also includes substances and substance mixtures of human origin. The centrifuge is only intended to be used for this purpose. A different use or application over and above this is deemed not in accordance with the specifications. The company Andreas Hettich GmbH & Co. KG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. Improper use and handling can result in dangers to life and limb of the user or third parties and impairments to the machine or to other material assets. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

3 Notes on safety

No claim under guarantee will be considered by the manufacturer unless the above instructions have been adhered to.

- **Before the initial operation of your centrifuge you should read and pay attention to the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device.**
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate.
 - However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge should be installed on a good, stable base.
- Before using the centrifuge absolutely check the rotor for firm placement.
- When the centrifuge is running, according to IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either run hot in the cold room for 30 minutes, or warm up for 3 hours in the warm room, before connecting to the mains.
- Load centrifuge rotor evenly. All positions on rotor must be filled.
- Do not fill centrifuge containers inside the centrifuge.
- Centrifuge containers must not be filled beyond the capacity specified by the manufacturer.
- Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2)
- Only the rotors and accessories approved by the manufacturer for this device may be used (see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories")."
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures. Without additional proceedings (like an additional bioseal between bucket and lid of bucket or angle rotor with a special bioseal between rotor and lid) a centrifuge is not a biosafety system in accordance to the regulation EN 61010-2-020. In the case of material belonging to risk group II (see the World Health Organisation's "Laboratory Biosafety Manual") they should employ a biosafety system. Under this system small drips and aerosols are prevented from escaping by a bioseal (packing ring) located between the hanger and the lid. Centrifuge containers with special screw caps, as obtainable through trade suppliers, can also be used for hazardous substances.

In the case of materials from the higher risk groups greater safety provision is required than the arrangements described above. In a biosafety system, centrifuge containers with special screw caps must be used.

- For further details of available biosafety systems see section "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.
- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Any rotors, hangers or accessories showing clear signs of corrosion or mechanical defects must not be used for centrifuging.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Hettich company are allowed to be utilised.
- The following safety regulations apply:
IEC 61010-1 and IEC 61010-2-020 as well as their national deviations..
- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of IEC stipulations,
 - prescribed tests to BGV A1, BGR 261 are carried out by an expert.

4 Symbol meanings

Symbol on the machine:

Attention, general hazard area.

Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!

Symbol in the operating instructions:

Attention, general hazard area.

This symbol refers to safety relevant warnings and indicates possibly dangerous situations.

The non-adherence to these warnings can lead to material damage and injury to personal.

This symbol refers to important circumstances.

Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).

Applies in the countries of the European Union, as well as in Norway and Switzerland.

5 Delivery checklist

The following items and accessories are delivered with the centrifuge:

- | | |
|---|--------------------------------------|
| 1 | Connecting cable |
| 2 | Fuses |
| 2 | Carbon brushes |
| 1 | Notes on moving the equipment safely |
| 1 | Operating instructions |

The model EBA 20 is supplied complete with angle-rotor 8x15 ml. For round-bottomed tubes no rubber inserts are necessary. Adapters for various types of blood collecting tubes can be offered on request.

6 Space requirement

- The necessary space requirement can be found under Dimensions in the Technical data chapter.
- The centrifuge must be set up in a suitable place, so that it is stable. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to IEC 61010-2-020.

When the centrifuge is running, according to IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

- Substances can be discharged through the ventilation opening on the rear side of the centrifuge. The device is to be placed in such a way that the air stream is not directed towards people.
- Do not place any object in front of the ventiduct.
 - Keep a ventilation area of 300 mm around the ventiduct.

7 Connection to the mains

- Check whether the supply voltage, supply frequency and on-site mains fuse agree with the specification on the nameplate. The nameplate is located on the back of the centrifuge.
- Make sure that the mains switch is in the “0” position.
- The centrifuge must be connected to a standard mains socket using the power supply cable provided.

8 Commissioning

- Check that the centrifuge has been properly set up and that the electrical connections are correct (see Connection to the Mains and Space Requirement).
- Switch the mains switch “ON”, switch position “I”.
After a short time the control panel will switch on and the last set parameters will appear in the displays.
- When the  symbol has lit up, the lid can be opened.
- Open the lid and remove the transport protection (see Transport Protection information sheet).

9 Opening and closing the lid

9.1 Opening the lid

- When the  symbol has lit up, the lid can be unlocked and opened.
- Unlock the lid by pressing the  key and swivel the lid upwards.
- The  symbol appears.

The lid can only be opened if the centrifuge is switched on and the rotor is stationary.
If this is not possible, see chapter “Emergency release”.

9.2 Closing the lid

- Close the lid by swinging the lid down and pressing lightly. If the symbol  lights up in the rotation display the lid has been correctly locked.

Do not bang the lid shut.


10 Arming the rotor

- Always fill the centrifuge containers outside the centrifuge.
- Ensure, by eye, that the containers are filled to a uniform level.
- Load positions located opposite one another equally.
- Ensure that the rotor is correctly positioned and secured.
- The weight of the permissible filling quantity is specified on each rotor. This weight may not be exceeded.

11 Control panel

The control panel is divided into four areas.

- Speed
- Rotation display
- Time
- Keys

11.1 Speed area

The speed area contains the speed display which indicates the preselected speed or the ACTUAL speed during operation.

The speed can be preselected or changed during operation with the ▲ and ▼ arrow keys. At the end of operation, the display flashes in seconds until the lid is opened or a key is activated.

11.2 Rotation area

The rotation area contains the rotation display, which illuminates and rotates anti-clockwise whilst the rotor is turning. When the  symbol lights up, the rotor is stationary and the lid can be opened.

When the lid is unlocked, the  symbol lights up.

11.3 Time area

This area contains the running time display and the ▲ and ▼ arrows, with which the switch-on time can be preselected or changed during operation.

When stationary, the preselected switch-on time is displayed, as well as the remaining run time during operation. The decimal point flashes each second during operation.

At the end of operation, the display flashes each second until the lid is opened or a key is activated.

11.4 Key area

Start key 	To start a run with the preselected parameters.
Stop key 	The STOP key allows a centrifuging run to be stopped at any time.
Pulse key 	Brief centrifuging whilst the key is pressed down. The run time is displayed in seconds in the run time display. After 60 seconds, the run time is displayed in minutes.
Lid key 	After standstill the lid can be unlocked by pressing the corresponding key.

12 Adjustable parameters

Speed	<ul style="list-style-type: none"> min. speed 500 rpm max. speed 6000 rpm in steps of 100 rpm Adjustable up to max. speed of the rotor used.
Time	<ul style="list-style-type: none"> min. preset time 1 min max. preset time 99 min in steps of 1 min or
"—" 	<ul style="list-style-type: none"> continuous operation or
Braking force	<ul style="list-style-type: none"> full braking effect low braking effect

13 Preselecting centrifuging parameters or changing them during operation

13.1 Speed

Each time the  and  arrow keys are activated, the speed is changed by 100 rpm. If an arrow key is kept pressed down for longer, the speed is changed at an increasing rate.

If the speed is changed during operation, this becomes effective immediately, i.e. the speed is adjusted.

13.2 - of denser substances

The rotors are designed to centrifuge substances up to a maximum mean homogenous density of 1.2 kg/dm³ when rotating at the stated speed.

Denser substances must be centrifuged at lower speed.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed (n}_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density}}} \times \text{Rated speed}$$

e.g.: RPM 4000, density 1.6 kg/dm³

$$n_{\text{red}} = \sqrt{\frac{1.2}{1.6}} \times 4000 = 3464 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

13.3 Run time

Each time the  and  arrow keys are activated, the run time is changed by 1 minute. If an arrow key is kept pressed down for longer, the run time is changed in minute steps, at an increasing rate.

If the run time is changed during operation, this becomes effective immediately, i.e. the remaining run time is adjusted.

If an arrow key is activated during operation, the run time or speed is adjusted accordingly.

13.3.1 Continuous operation

The centrifuge can be operated in continuous operation. Continuous operation is set by pressing the  arrow key until "—" appears in the run time display. Operation is only stopped by activating the  key.


13.3.2 Pulsed operation


For brief centrifuging. The rotor turns at the preselected speed as long as the  key is pressed down.


13.4 Brake adjustment

With this centrifuge, the braking effect can be set to normal or low.

This can be set before a run, as follows:

1. Switch off the mains switch
2. Simultaneously press the pulse key  and the arrow key  in the speed area and switch on the mains switch. Release the keys again.
3. If required press the arrow button  in the speed field until the depicted display appears. In the speed field the machine version as set by the factory and in the time field the brake setting.

4. Set the required value, "1" or "0", in the time area with the  and  arrow keys.
5. Confirm the set value with the  key. See Rotors and Accessories in the Braking Times chapter.

13.5 Start centrifuging run

 The centrifugation run can be interrupted at any time by pressing the key . The time and speed can be changed during the centrifugation run, with the keys  . If the key  or  is kept pressed, the value changes with increasing speed. After a centrifugation run, the display flashes until the cover is opened or a key is pressed. If the symbol "—" (lid closed) and "L" (lid open) flashes alternately in the rotation indicator , operation of the centrifuge can only be continued after opening the lid.

After setting the centrifuging parameters, start the centrifuge by pressing the  key.

13.6 End centrifuging run

A centrifuging run is ended by pressing the  key or after expiry of the set run time. Once the rotor has come to a stop the lid lock is briefly triggered and the display flashes. The lid can now be opened.

14 Calculating rotational speed RPM and relative centrifugal force RCF

These values are calculated using the formulas below:

$$RCF = \left(\frac{RPM}{1000} \right)^2 \times r \times 1,118 \qquad RPM = \sqrt{\frac{RCF}{r \times 1,118}} \times 1000$$

RCF = Relative centrifugal force
 RPM = Rotational speed (revolutions per minute)
 r = Radius in mm = Distance from the centre of the axle to the floor of the centrifuge container. For further details on radius see "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories" section.


15 Emergency release

The lid cannot be opened during power failure. An emergency release has to be executed by hand.

For emergency release disconnect the centrifuge from the mains.
Open the lid only during rotor standstill.
Only the plastic release pin provided may be used for emergency release.

1. Insert the release pin horizontally into the hole located in the middle of the front panel, see drawing.
2. Push the release pin to the right, then up to a 45° position. Pull the lid upwards at the same time to unlock, see drawing.
3. Remove the release pin from the hole again.

16 Maintenance and servicing

Pull the mains plug before cleaning.
Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Cleaning agents and disinfectants which lie in the pH range 5 – 8 are to be utilised. Alkaline cleaning agents with a pH value > 8 are to be avoided.
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

16.1 Centrifuge

- Regularly clean the centrifuge housing and the centrifugal chamber; and if necessary, clean with soap or a mild cleaning agent and water. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately. For surface disinfection we recommend Bacillol Plus from the company Bode Chemicals in Hamburg. After use the disinfection medium Bacillol Plus must be thoroughly removed with water.
- Lightly grease the rubber seal of the centrifugal chamber after every cleaning.

16.2 Rotors and Attachments

- Rotors and accessory parts must be regularly cleaned with soap or a mild cleaning agent and water in order to prevent corrosion and changes of material. Cleaning is recommended at least once a week, even better after every usage.
- If the rotor or accessory parts are contaminated by pathogenic or radioactive material, a suitable cleaning has to be executed. For disinfection we recommend Helipur H plus N from the company B. Braun Melsungen. For the removal of radioactive material we recommend deco neutracon from the company Decon Laboratories Limited.
- The rotors and accessory parts must be dried immediately after cleaning.
- The rotors and accessory parts are to be checked on a monthly basis for corrosion damage.

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

16.3 Autoclaving

The rotor and accessories must not be autoclaved.

16.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splitters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

17 Faults

If the fault cannot be rectified according to the faults table, Hettich customer services must be informed.

Please state the type of centrifuge and the factory serial number. Both values are visible on the centrifuge type plate.


MAINS RESET: - Mains switch OFF for longer than 10 secs.
- Mains switch ON.

Fault	Display	Cause of fault	Remedy
No display	---	No voltage Triggering of the excess current cut-out.	- Check distribution voltage. - Check fuses. - Mains switch ON.
Tacho error	- 1 -	Failure of speed impulses during operation, Unbraked runout.	- Wait 100 sec. delay time (lid must be closed). - Perform mains reset.
System reset	- 2 -	Failure of power supply during a run Unbraked runout.	- When stationary, open lid and press (START) key.
Kommunikation	- 4 -	Fault in control unit or power unit. Unbraked runout.	- Perform mains reset when stationary.
Overload	- 5 -	Motor or motor control defective.	- Perform mains reset when stationary.
Overspeed	- 7 -	Fault in the supply board Unbraked runout	- Perform mains reset when stationary
Lid error	- d -	Unbraked rundown, after standstill lid release.	- Perform mains reset after standstill.

18 Change mains input fuse

Switch off the mains switch and separate the centrifuge from the mains!

The fuse holder (A) with the mains input fuses is located next to the mains switch.

- Remove the connecting cable from the machine plug socket.
- Press the snap-fit (B) against the fuse holder (A) and remove.
- Exchange defective mains input fuses.

Only use fuses with the rating defined for the type. See the following table.

- Reinsert the fuse holder until the snap-fit clicks shut.
- Reconnect the centrifuge to the mains supply.

Model	Type	Fuse	Order no.
EBA 20	2002	T 1,6 AH/250V	E891
EBA20	2002-01	T 3,15 AH/250V	E997

19 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

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All rights reserved. No part of this publication may be reproduced without the written prior permission of the copyright owner.